# Airsonic Mini-Sandblaster

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# Airsonic Mini-Sandblaster

## Mini-Sandblaster for intra and extra oral use

The **AIRSONIC MINI-SANDBLSTER** greatly enhances the adherence of cements by enlarging the size of the accordingly treated surfaces.

### The **AIRSONIC MINI-SANDBLSTER**

- optimizes the retention to and the adherence of a great variety of materials
- effectively reduces the number of tries for a solid adherence
- is easy to use and to install
- is the ideal aid for the adhesive technique in the dental office and laboratory because of its substantially improved retention / adherence to most surfaces such as:
  - precious metal as well as non-precious metal
  - porcelain
  - composites
  - resins
  - hybrid micro fillings

### Installation of the **AIRSONIC MINI-SANDBLSTER?**

The sandblaster requires compressed air (4 - 8 bar). In case that dry air compressor is not available, a water trap filter must be installed using a tee fitting and a connecting tube. (The technical service of your dental dealer will assist you if necessary):

1. Cut the air supply line and install the tee fitting.
2. Connect the female part of the fitting, which contains the automatic shut off, to the tee fitting.
3. Turbine connection kits for Kavo, Sirona-Titane or Roroquick (W&H) optionally available. The kit connects the Airsonic unit to the existing airtube by means of a T.connecting piece.
4. Adapters for the laboratory connections and EMS box (old Type) are optionally available also.

### Operations of the **AIRSONIC MINI-SANDBLSTER?**

When the air supply line has been connected to the unit as described above the sandblaster is ready for use. For abrading please use **AIRSONIC Alu-Oxide Powder** with 50 micron (order no. 605 084) or 90 micron (no. 605 086). The powder should be absolutely dry and fine-grain. Fill the transparent jar on the handpiece 3/4 full.

The sandblaster is held like a pen using the thumb to activate the finger button control. The distance between nozzle and the surface should be 2 - 10 mm. Light pressure on the control button will activate the sandblasting. The strenght of the abrasive jet stream is increased by a firmer pressure on the control button. The nozzle can be rotated 360° from side to side to obtain the optimal angle of the nozzle to the object. If your sandblaster should become clogged due to moist powder, press the button several times for short intervals. Push the nozzle with your finger as much sideways as possible to loosen clogging particles.

### Warranty:

The producer ensures that the product is free of defects and that it shows all guaranteed features and/or all characteristics agreed upon. The warranty starts with the day of delivery of the product and is valid for twelve months. Excluded from the warranty are faults, careless use of the product, connection to unsuitable power supply, operation with the wrong voltage, acts of God, such as fire, lightning, humidity, etc. In case of justified complaints the producer has the right to repair, to deliver spare parts or to decrease the sales price. Cancellation of the sales contract can only be claimed by the customer after two repairs have been done without success.