


THE NEW "Super" Quick-Tip

More Precise
More Consistent
More Resistant



1/12 of a drop each
and every time



"Super" Quick-Tip

Other Micro Applicators

- For application of minute amounts of liquids and pastes: cavity lining, etching, bonding, fluoride varnishing, hemostat and disclosing solutions.
- Non-linting, non-absorbing tips consistently release precise amounts of liquids.
- Flexible tapered ends maintain desired angle for application in difficult areas.

The product of an extensive research program, "Super" Quick-Tips retain all the benefits of the regular Quick-Tips, plus are completely resistant to high concentrations of phosphoric acid and acetone. Testing by an independent laboratory shows that competitive fiber micro-applicators may "melt" in these solutions and leave a residue on the surface of the tooth. "Super" Quick-Tips are recommended as the only micro applicator tip for use with all etching and bonding materials. Super Quick-Tips now come in a hygienic white color available in dispenser packages of 100 or 250 tips.



To place an order, call your preferred dental dealer.
For a dealer near you call:
(800) 328-2335
Fax: **(800) 573-9392**



HAGER
WORLDWIDE

12890 Automobile Blvd., Suite A
Clearwater, FL 33762