

**NOW
AVAILABLE**

BullFrog Invincible



TITANIUM


High Volume Evacuation Handpieces & Saliva Ejectors

Completely Resistant to Steam Autoclave, Dry Heat, Chemiclave, Liquid Sterilants, Ultra-Sonic Solutions, Enzymatic Solutions, Holding Solutions and *Sodium Hypochlorite* (Bleach) .

Available Lengths & Styles:

Ten (10) Year Warranty


Short



Standard Invincible HVE Handpiece

Titanium Body with
1/2"/1.2cm Swivel Tubing Adapter (Aluminum)
10 Yr. Warranty
Item# B11105T


Medium



Uni-Vac Invincible HVE Handpiece

Titanium Body with
1/2"/1.2cm Swivel Tubing Adapter (Aluminum)
10 Yr. Warranty
Item# B11150T


Long



Extended Invincible HVE Handpiece

Titanium Body
with 1/2"/1.2cm Swivel Tubing Adapter (Aluminum)
10 Yr. Warranty
Item# B11106T

Saliva Ejector



Solid Invincible Saliva Ejector

Titanium Body with
1/4"/.6cm Swivel Tubing Adapter (Aluminum)
10 Yr. Warranty
Item# B22257T

Bull Frog Invincible units are available singly with or w/o 1/2", 5/8", or Cox Titanium Adaptors or in a Combo Pak with a set of three at a 15% savings.

See HVE Tips on Reverse Side



HÄGER

WORLDWIDE

12890 Automobile Blvd. Clearwater, FL 33762
Tel: (800) 328-2335 (813) 573-1422
Fax: (800) 573-9392 (813) 541-3496